



Product: Anti-Bovine Osteocalcin
Catalog #: ABOC-5021
Lot Number: N1122A
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 9.5 mg/ml
Specific Activity: NA
Recommended Storage: -20°C
Buffer: 50% Glycerol/H₂O (v/v)
Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Mouse monoclonal IgG(1); antigen recognized: human and bovine bone osteocalcin; RIA(+), Blot(+), ELISA(+), Immunohistochemistry(+); calcium dependent.

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!