



Product: Anti-Bovine Factor V
Catalog #: ABV-5106
Lot Number: J0731
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 9.4 mg/ml
Specific Activity: NA
Recommended Storage: -20°C
Buffer: 50% Glycerol/H₂O (v/v)

Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

-Mouse monoclonal IgG(1); Antigen recognized: heavy chain region of bovine factors V and Va; RIA/EIA (+), Blot (+), Immunofluorescence (+).