



Product:	Anti-Bovine Factor X
Catalog #:	ABX-5051
Lot Number:	R0429
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	5.2 mg/ml
Specific Activity:	NA
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H <sub>2</sub> O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

**Product/Caution Notes:**

Technical and use information is available at [www.haemtech.com/FAQ.htm](http://www.haemtech.com/FAQ.htm)

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

-Mouse monoclonal IgG(1); Antigen recognized: heavy chain regions of factors X and Xa (reactive toward human, bovine, rabbit, sheep, porcine and canine factor X); Useful for purification.  $K_d=9 \times 10^{-11}$ M, RIA(+), Blot (+). Inhibitory, partial calcium dependence.