



Product:	Anti-Human Factor IX
Catalog #:	AHIX-5041
Lot Number:	T0830
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	7.1 mg/ml
Specific Activity:	N/A
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H ₂ O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

*Affinity Purified.

Mouse monoclonal IgG(1); antigen recognized: Human factor IX, human factor IXa, heavy chain of human factors IX and IXa; RIA(+), Blot(+), ELISA(+), Immunohistochemistry(+).

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!