



Product: Anti-Human Protein C  
Catalog #: AHPC-5071  
Lot Number: X1117  
Extinction Coefficient:  $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$   
Molecular Weight: 150000  
Concentration: 8.8 mg/mL  
Specific Activity: ELISA (+)  
Recommended Storage: -20 °C  
Buffer: 50% Glycerol/H<sub>2</sub>O (v/v)  
Amount Shipped: \_\_\_\_\_  
Date Shipped: \_\_\_\_\_  
Shipped By: \_\_\_\_\_  
Expiration: \_\_\_\_\_

**Product/Caution Notes:**

Technical and use information is available at [www.haemtech.com/FAQ.htm](http://www.haemtech.com/FAQ.htm)

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Mouse monoclonal IgG(1); Antigen recognized: Human protein C, and aPC. ELISA(+), RIA (+), Blot (+); Useful for purification; Partial calcium dependence.