



Product: Anti-Human Protein S
Catalog #: AHPS-5091
Lot Number: Y0220
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 4.0 mg/mL
Specific Activity: N/A
Recommended Storage: -20 °C
Buffer: 50% Glycerol/H₂O (v/v)

Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Mouse monoclonal IgG(2b); Antigen recognized: human protein S and Protein S-C4bp complex; Useful for purification. RIA(+), Blot (+) non-reduced, calcium dependent, ELISA (+).

-ELISA(+).