



Product: Anti-Human Prothrombin
Catalog #: AHP-5013
Lot Number: U1213
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 10.2 mg/mL
Specific Activity: N/A
Recommended Storage: -20 °C
Buffer: 50% Glycerol/H₂O (v/v)

Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Mouse monoclonal IgG(2a); antigen recognized: Human prothrombin, prethrombin-1, fragment 2, meizothrombin; ELISA(+), Blot(+); inhibits clotting and prothrombin activation.

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!