



Product:	Anti-Human TAFI
Catalog #:	AHTAFI-5024
Lot Number:	N0531
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	10 mg/ml
Specific Activity:	NA
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H ₂ O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Mouse monoclonal IgG(1); Antigen recognized: Human TAFI and activated TAFI; Blot (+), ELISA (+); Inhibits TAFI activation; Inhibits activated TAFI.