



Product: Anti-Human TAFI
Catalog #: AHTAFI-5065
Lot Number: N0425
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 10.8 mg/mL
Specific Activity: N/A
Recommended Storage: -20°C
Buffer: 50% Glycerol/H₂O (v/v)
Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Mouse monoclonal IgG(3); Antigen recognized: Human TAFI, activated TAFI. Blot (-), ELISA (+); Does not inhibit TAFI activation; Inhibits activated TAFI.