



Product: Anti-Human Tissue Factor
Catalog #: AHTF-5264
Lot Number: T0503
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 4.7 mg/ml
Specific Activity: N/A
Recommended Storage: -20°C
Buffer: 50% Glycerol/H₂O (v/v)

Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

Mouse monoclonal IgG; antigen recognized: Human tissue factor; Blot(+), ELISA(+).