



Product: Anti-Human Thrombin
Catalog #: AHT-5020
Lot Number: U0825
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 7.8 mg/ml
Specific Activity: N/A
Recommended Storage: -20°C
Buffer: 50% Glycerol/H₂O (v/v)
Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Mouse monoclonal IgG(1); Antigens recognized: human thrombin (Kd = 1.4 X 0.0000001), thrombin-ATIII complex (Kd = 1.5 X 0.0000001), thrombin-FPRck; ELISA (+), Blot (-).