



Product: Anti-Human Factor V
Catalog #: AHV-5101
Lot Number: W0725
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 13.9 mg/mL
Specific Activity: N/A
Recommended Storage: -20 °C
Buffer: 50% Glycerol/H₂O (v/v)

Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Mouse monoclonal IgG(1); Antigen recognized: light chain of human factor V, human factor V, human factor Va, bovine factor Va; K_d=3X10⁽⁻⁹⁾M, RIA(+), Blot(-), inhibitory.

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!