



Product: Anti-Human Factor V
Catalog #: AHV-5102
Lot Number: G0517
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 15.4 mg/ml
Specific Activity: NA
Recommended Storage: -20°C
Buffer: 50% Glycerol/H₂O (v/v)

Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Mouse monoclonal IgG; Antigen recognized: 120,000 activation peptide of human factor V, not bovine factor V or Va; $K_d=4 \times 10^{-9}\text{M}$, RIA(+), Blot(+).

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!