



Product: Anti-Human Factor V
Catalog #: AHV-5110
Lot Number: L0930
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 9 mg/ml
Specific Activity: NA
Recommended Storage: -20 °C
Buffer: 50% Glycerol/H₂O (v/v)

Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Mouse monoclonal IgG(1); antigen recognized: 120,000 dalton activation peptide of human factor V; RIA(+), Blot(+); useful for the purification of the activation peptide.

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!