



Product:	Anti-Human Factor X
Catalog #:	AHX-5050
Lot Number:	W1214
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	11.6 mg/mL
Specific Activity:	N/A
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H ₂ O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

-Mouse monoclonal IgG(1); Antigen recognized: human factor X, heavy chain regions of human factors X and Xa; Useful for purification; Inhibits PT, prothrombinase and partially aPTT, but not amidase activity.