



Product:	Rat Anti-Mouse Factor IXa
Catalog #:	AMIXA-9042
Lot Number:	X0116
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	7.5 mg/mL
Specific Activity:	N/A
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H <sub>2</sub> O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

**Product/Caution Notes:**

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

Technical and use information is available at [www.haemtech.com/FAQ.htm](http://www.haemtech.com/FAQ.htm)

ELISA (+) for mouse IX and IXa. Western Blot (+) for mouse IX (Non-reduced and weak reduced) and Ixa (non-reduced only). Non-inhibitory in aPTT Assay.