



Product: Rat Anti-Mouse Prothrombin  
Catalog #: AMP-9013  
Lot Number: W0326  
Extinction Coefficient:  $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$   
Molecular Weight: 140000  
Concentration: 4.0 mg/mL  
Specific Activity: NA  
Recommended Storage:  $-20^{\circ}\text{C}$   
Buffer: 50% Glycerol/H<sub>2</sub>O (v/v)

Amount Shipped: \_\_\_\_\_  
Date Shipped: \_\_\_\_\_  
Shipped By: \_\_\_\_\_  
Expiration: \_\_\_\_\_

**Product/Caution Notes:**

ELISA (+) for mouse prothrombin. Blot (+) for non-reduced mouse prothrombin.

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

Technical and use information is available at [www.haemtech.com/FAQ.htm](http://www.haemtech.com/FAQ.htm)