



Product:	Rat Anti-Mouse Factor VIII
Catalog #:	AMVIII-9035
Lot Number:	W0612
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 1.4$
Molecular Weight:	150000
Concentration:	8.7 mg/mL
Specific Activity:	N/A
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H ₂ O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

Technical and use information is available at www.haemtech.com/FAQ.htm

Non-inhibitory in aPTT; ELISA (+) for rmFVIII, no data formFVIIIa; Blot (+) rmFVIII (non-reduced only), no data for mFVIIIa, hFVIII, hFVIIIa.