



Product:	Rat Anti-Mouse Factor VII
Catalog #:	AMVII-9031
Lot Number:	W0920
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	8.7 mg/mL
Specific Activity:	N/A
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H ₂ O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

Technical and use information is available at www.haemtech.com/FAQ.htm

ELISA (+) for recombinant mFVII/FVIIa; native mFVII/FVIIa not determined. Western (+) for recombinant mFVII/FVIIa & native mFVII (non-reduced only); native mFVIIa not determined; does not cross-react with hFVII/FVIIa. Non-inhibitory in a PT-based clotting assay.