



Product:	Rat Anti-Mouse Factor X
Catalog #:	AMX-9050
Lot Number:	W0302
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	8.0 mg/mL
Specific Activity:	na
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H ₂ O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

Technical and use information is available at www.haemtech.com/FAQ.htm

ELISA (+) for mouse X and Xa. Blot (+) for mouse X and Xa heavy chain (reduced/non-reduced), and cross reacts with human X and Xa heavy chain (non-reduced - strong; reduced - weak). Does not inhibit PT assay (mouse plasma).