



Product: Rat Anti-Mouse Factor X
Catalog #: AMX-9051
Lot Number: W0302
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 8.0 mg/mL
Specific Activity: na
Recommended Storage: -20 °C
Buffer: 50% Glycerol/H₂O (v/v)
Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

Technical and use information is available at www.haemtech.com/FAQ.htm

ELISA (+) for mouse X, but not Xa. Blot (+) for mouse X heavy chain (reduced - strong; non-reduced - weak), and does NOT cross react with human X or Xa. Does not inhibit PT assay (mouse plasma).