



Product:	Anti-Human Osteonectin
Catalog #:	AON-5031
Lot Number:	W1214
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	10.8 mg/mL
Specific Activity:	N/A
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H ₂ O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

Antibody-Antigen interaction is calcium dependent.

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Mouse monoclonal IgG(1); Antigen recognized: Human, rat and bovine bone osteonectin, human platelet osteonectin, SPARC produced by various cells grown in culture; RIA (+), Blot (+), ELISA (+); Immunohistochemistry (+); Kd(app) = 5 X 0.00000001 M.