



Product:	DAPA*
Catalog #:	DAPA
Lot Number:	X0512-1MG
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = \text{N/A}$
Molecular Weight:	539
Concentration:	9.4 mg/mL (17.4 mM)
Specific Activity:	N/A
Recommended Storage:	at or below -70 °C
Buffer:	H ₂ O
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

This product is light sensitive. Protect from light!

*Dansylarginine N-(3-ethyl-1,5-pentanediy)amide.