



Product: _____ HUMAN FACTOR XIII A
Catalog #: _____ HCXIIIA-0165
Lot Number: _____ X0506-0.05MG
Concentration: _____ 2.8 mg/mL
Specific Activity: _____ 345 U/mg*
Extinction Coefficient: _____ $E_{1\text{cm}, 280\text{nm}}^{1\%} = 13.8$
Molecular Weight: _____ 312000
Buffer: _____ 50% Glycerol / 500 uM EDTA (vol/vol)
Recommended Storage: _____ -20 ° C
Amount Shipped: _____
Date Shipped: _____
Shipped by: _____
Expiration: _____

PRODUCT/CAUTION NOTES:

Products contained within this shipment are of human blood/plasma origin. Although the starting material was tested prior to initiation of the manufacturing process, and was found negative or nonreactive for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II, extreme caution should be used when handling this material as there is a margin of error in all tests.

This product is to be used for IN VITRO research purposes only, and is not to be used for clinical, diagnostic or human use!

*1 unit is equal to the Factor XIII activity in 1 milliliter of normal plasma.