



Product: Anti-Human Protein C
Catalog #: AHPC-5071
Lot Number: X1117
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight: 150000
Concentration: 8.8 mg/mL
Specific Activity: ELISA (+)
Recommended Storage: -20 °C
Buffer: 50% Glycerol/H₂O (v/v)
Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Mouse monoclonal IgG(1); Antigen recognized: Human protein C, and aPC. ELISA(+), RIA (+), Blot (+); Useful for purification; Partial calcium dependence.