



Product:	Anti-Human Factor V
Catalog #:	AHV-5101
Lot Number:	W0725
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	13.9 mg/mL
Specific Activity:	N/A
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H <sub>2</sub> O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

**Product/Caution Notes:**

Mouse monoclonal IgG(1); Antigen recognized: light chain of human factor V, human factor V, human factor Va, bovine factor Va; Kd=3X10(-9)M, RIA(+), Blot(-), inhibitory.

Technical and use information is available at [www.haemtech.com/FAQ.htm](http://www.haemtech.com/FAQ.htm)

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!