



Product: Anti-Human Factor XI  
Catalog #: AHXI-5061  
Lot Number: X0116  
Extinction Coefficient:  $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$   
Molecular Weight: 150000  
Concentration: 12.4 mg/mL  
Specific Activity: N/A  
Recommended Storage: -20 °C  
Buffer: 50% Glycerol/H<sub>2</sub>O (v/v)

Amount Shipped: \_\_\_\_\_  
Date Shipped: \_\_\_\_\_  
Shipped By: \_\_\_\_\_  
Expiration: \_\_\_\_\_

**Product/Caution Notes:**

Technical and use information is available at [www.haemtech.com/FAQ.htm](http://www.haemtech.com/FAQ.htm)

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Mouse monoclonal IgG; Antigen recognized: Human factor XI and XIA. Useful for purification; Blot (+) nonreduced only; RIA (+). Inhibitory in clotting assay.