



Product: Rat Anti-Mouse Protein C
Catalog #: AMPC-9071
Lot Number: U0816
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 1.4$
Molecular Weight: 150000
Concentration: 9.0 mg/mL
Specific Activity: N/A
Recommended Storage: -20 °C
Buffer: 50% Glycerol/H₂O (v/v)

Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

-SPECIFICITY: ELISA (+), PC only (not aPC). Western Blot (+), PC only (not aPC), does not cross-react with human PC/aPC. Does not inhibit aPC. Inhibits PC activation.

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!