



Product:	Rat Anti-Mouse Protein C
Catalog #:	AMPC-9072
Lot Number:	U0823
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14$
Molecular Weight:	150000
Concentration:	7.7 mg/mL
Specific Activity:	N/A
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H <sub>2</sub> O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

**Product/Caution Notes:**

Technical and use information is available at [www.haemtech.com/FAQ.htm](http://www.haemtech.com/FAQ.htm)

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

-SPECIFICITY: ELISA (+), PC and aPC. Western Blot (+), PC only (not aPC), does not cross-react with human PC/aPC. Does not inhibit aPC. Weak inhibition of PC activation.