



Product: Rat Anti-Mouse Factor X  
Catalog #: AMX-9051  
Lot Number: W0302  
Extinction Coefficient:  $E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$   
Molecular Weight: 150000  
Concentration: 8.0 mg/mL  
Specific Activity: na  
Recommended Storage:  $-20^{\circ}\text{C}$   
Buffer: 50% Glycerol/H<sub>2</sub>O (v/v)

Amount Shipped: \_\_\_\_\_  
Date Shipped: \_\_\_\_\_  
Shipped By: \_\_\_\_\_  
Expiration: \_\_\_\_\_

**Product/Caution Notes:**

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

Technical and use information is available at [www.haemtech.com/FAQ.htm](http://www.haemtech.com/FAQ.htm)

ELISA (+) for mouse X, but not Xa. Blot (+) for mouse X heavy chain (reduced - strong; non-reduced - weak), and does NOT cross react with human X or Xa. Does not inhibit PT assay (mouse plasma).