



Product:	BETA-2-GLYCOPROTEIN I
Catalog #:	B2GI-0001
Lot Number:	Y0819-0.1MG
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 10.0$
Molecular Weight:	54200
Concentration:	Lyophilized*
Specific Activity:	4.51 Plasma Equalency Units per milligram
Recommended Storage:	4 °C
Buffer:	0.2 M Glycine, 0.15 M NaCl, pH 7.4
Amount Shipped:	_____
Date Shipped:	_____
Shipped By:	_____
Expiration:	_____

**Product/Caution Notes:**

Products contained within this shipment are of human blood/plasma origin. Although the starting material was tested prior to initiation of the manufacturing process, and was found negative or nonreactive for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II, extreme caution should be used when handling this material as there is a margin of error in all tests.

This product is to be used for IN VITRO research purposes only, and is not to be used for clinical, diagnostic or human use!

\*Add 1 milliliter of water for a final concentration of 0.1 milligram/milliliter.

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

-ELISA (+).