



Product: HUMAN ALPHA-2-ANTIPLASMIN
Catalog #: HA2AP-0230
Lot Number: Y0120-0.1MG
Concentration: 8.5 mg/mL
Specific Activity: 0.83 mol HCPM/ mol a-2*
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 7.03$
Molecular Weight: 70000
Buffer: 50mM Potassium Phosphate, 7.5mM Potassium Chloride, 75uM EDTA, pH 7.4
Recommended Storage: at or below -70 ° C
Amount Shipped: _____
Date Shipped: _____
Shipped by: _____
Expiration: _____

PRODUCT/CAUTION NOTES:

*Specific activity indicates the amount of human alpha-2-antiplasmin (on a molar basis) that inhibits 1 mole of human plasmin.

Products contained within this shipment are of human blood/plasma origin. Although the starting material was tested prior to initiation of the manufacturing process, and was found negative or nonreactive for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II, extreme caution should be used when handling this material as there is a margin of error in all tests.

This product is to be used for IN VITRO research purposes only, and is not to be used for clinical, diagnostic or human use!