



Product: Human Factor XIII
Catalog #: HCXIII-0160
Lot Number: HH0709-0.1MG
Extinction Coefficient: $E_{1\text{cm}, 280\text{nm}}^{1\%} = 13.8$
Molecular Weight: 320000
Concentration: 6.9 mg/mL
Specific Activity: 64 U/mg*
Recommended Storage: -20 °C
Buffer: 50% Glycerol/500 μM EDTA (vol/vol)

Amount Shipped: _____
Date Shipped: _____
Shipped By: _____
Expiration: _____

Product/Caution Notes:

Products contained within this shipment are of human blood/plasma origin. Although the starting material was tested prior to initiation of the manufacturing process, and was found negative or nonreactive for anti-HIV-1/2, HIV-1 RNA, HBsAg, STS, anti-HCV, HCV RNA, anti-HBcore and anti-HTLV I & II, and Zika (by investigational NAT) extreme caution should be used when handling this material as there is a margin of error in all tests, and testing does not always disclose pathogens or other harmful substances contained in blood/plasma products.

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION! In case of contact or ingestion consult a physician.

*1 unit is equal to the Factor XIII activity in 1 milliliter of normal plasma.