



|                         |                                             |
|-------------------------|---------------------------------------------|
| Product:                | Human Heparin Cofactor II                   |
| Catalog #:              | HCII-0190                                   |
| Lot Number:             | Y0302                                       |
| Extinction Coefficient: | $E_{1\text{cm}, 280\text{nm}}^{1\%} = 5.93$ |
| Molecular Weight:       | 65600                                       |
| Concentration:          | 7.4 mg/mL                                   |
| Specific Activity:      | 1220 U/mg                                   |
| Recommended Storage:    | -20 °C                                      |
| Buffer:                 | 50% Glycerol/ H2O (vol/vol)                 |
| Amount Shipped:         |                                             |
| Date Shipped:           |                                             |
| Shipped By:             |                                             |
| Expiration:             |                                             |

**Product/Caution Notes:**

Products contained within this shipment are of human blood/plasma origin. Although the starting material was tested prior to initiation of the manufacturing process, and was found negative or nonreactive for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II, extreme caution should be used when handling this material as there is a margin of error in all tests.

This product is to be used for IN VITRO research purposes only, and is not to be used for clinical, diagnostic or human use!

This product is greater than or equal to 95% Pure by SDS-PAGE

\*One Unit will inhibit one (1) NIH Unit of Thrombin.