



Product:	Human Fibronectin
Catalog #:	HCFN-0170
Lot Number:	X0606-2MG
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	550000
Concentration:	Lyophilized at 2.0 mg/mL+
Specific Activity:	N/A
Recommended Storage:	4 °C
Buffer:	20 mM Hepes, 150 mM NaCl, pH 7.4
Amount Shipped:	_____
Date Shipped:	_____
Shipped By:	_____
Expiration:	_____

**Product/Caution Notes:**

Products contained within this shipment are of human blood/plasma origin. Although the starting material was tested prior to initiation of the manufacturing process, and was found negative or nonreactive for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II, extreme caution should be used when handling this material as there is a margin of error in all tests, and testing does not always disclose pathogens or other harmful substances contained in blood/plasma products.

This product is to be used solely for IN VITRO research purposes only, and is not to be used for clinical, diagnostic or human use! In case of contact or ingestion consult a physician.

+Reconstitute with 1.02 milliliters of water to achieve 2.0 mg/mL.