



Product:	Human Protein C
Catalog #:	HCPC-0070
Lot Number:	Y0312-0.1MG
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.5$
Molecular Weight:	62000
Concentration:	7.6 mg/mL
Specific Activity:	<1% APC Activity*
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/H <sub>2</sub> O (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

**Product/Caution Notes:**

Products contained within this shipment are of human blood/plasma origin. Although the starting material was tested prior to initiation of the manufacturing process, and was found negative or nonreactive for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II, extreme caution should be used when handling this material as there is a margin of error in all tests.

This product is to be used for IN VITRO research purposes only, and is not to be used for clinical, diagnostic or human use!

\*Determined by chromogenic assay.