



Product: "ANGIOSTATIN"  
Catalog #: ANG-01  
Lot Number: U0502-0.5MG  
Concentration: 3.8 mg/mL  
Specific Activity: N/A  
Extinction Coefficient:  $E_{1\text{cm}, 280\text{nm}}^{1\%} = 17.4$   
Molecular Weight: 50000  
Buffer: 20 mM HEPES, 150 mM Sodium Chloride, pH 7.4  
Recommended Storage: -80 ° C  
Amount Shipped: \_\_\_\_\_  
Date Shipped: \_\_\_\_\_  
Shipped by: \_\_\_\_\_  
Expiration: \_\_\_\_\_

**PRODUCT/CAUTION NOTES:**

Products contained within this shipment are of human blood/plasma origin. Although the starting material was tested prior to initiation of the manufacturing process, and was found negative or nonreactive for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II, extreme caution should be used when handling this material as there is a margin of error in all tests.

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use!