



Product:	Horse Anti-Human Factor V
Catalog #:	PAHFV-H
Lot Number:	X1117
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	10.5 mg/mL
Specific Activity:	ELISA+
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/ H2O (vol/vol)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Specificity: Factor V, factor Va and their activation products; ELISA (+), Western blot (+), RIEP ().