



Product:	Sheep Anti-Human Thrombin
Catalog #:	PAHT-S
Lot Number:	M0225B
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	18 mg/ml
Specific Activity:	NA
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/10 mM Hepes, 50 mM NaCl, pH 7.4
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

Technical and use information is available at www.haemtech.com/FAQ.htm

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Specificity: Purified thrombin and prothrombin; ELISA (+), Western blot (), RIEP (); Useful for purification.