



Product:	Sheep Anti-Murine Factor X, Immunoabsorbed
Catalog #:	PAMFX-SIA
Lot Number:	P1106
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 14.0$
Molecular Weight:	150000
Concentration:	15.3 mg/ml
Recommended Storage:	-20 °C
Buffer:	50% Glycerol/20 mM Tris, 150 mM NaCl, pH 7.4 (v/v)
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

**Product/Caution Notes:**

Technical and use information is available at [www.haemtech.com/FAQ.htm](http://www.haemtech.com/FAQ.htm)

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

- Specificity: Mouse factor X; ELISA (+), Western blot (+), RIEP ( ); No cross reactivity with Human factor X (immunoabsorbed).