



Product:	Recombinant Tissue Factor
Catalog #:	RTF-0300
Lot Number:	Z0224
Extinction Coefficient:	$E_{1\text{cm}, 280\text{nm}}^{1\%} = 12.6$
Molecular Weight:	35000
Concentration:	1.5 mg/mL
Specific Activity:	3340 U/mg*
Recommended Storage:	at or below -70 °C
Buffer:	20 mM Tris, 150 mM NaCl, 10 mM CHAPS, pH 8.0
Amount Shipped:	
Date Shipped:	
Shipped By:	
Expiration:	

Product/Caution Notes:

This product is to be used for IN VITRO research purposes only, and is not intended for clinical or diagnostic use! NOT FOR HUMAN OR ANIMAL CONSUMPTION!

*Specific Activity determined by flourogenic assay.